


Table of Contents

	Page
Diagram: Instrument Classification Level 1 - 2	2
Instrument Classification Overview Level 1 - 3	3
Instrument Grouping	6
- Clinical Chemistry Instruments	6
- Immunochemistry Instruments	14
- Haematology Instruments	17
- Microbiology Instruments	22
- Other Clinical Instruments	23
- Data Management Systems	23
Alphabetical Order by Company	24
Alphabetical Order by Instrument Name	36

Diagram Instrument Classification Level 1 - 2



Level 3 Subgroups ==> see following pages

EDMA CODE	Company	Instrument-Name	Page
		<u>Instrument Classification Overview</u>	
		<u>Level 1 - 3</u>	
2		INSTRUMENTS	
21		CLINICAL CHEMISTRY INSTRUMENTS	6
21 01		Chemistry Analyser	
21 01 01		Manual Photometer and Spectrophotometer	
21 01 11		Dedicated Glucose Systems	7
21 01 12		Dedicated Glycated Haemoglobin /HbA1 Systems	
21 01 19		Other Dedicated Single-Parameter CC Systems	
21 01 29		Dedicated Multi-Parameter CC Systems	
21 01 31		Small Manual and Semi-automated CC Systems	
21 01 32		Semi-automated CC Analyser	
21 01 41		Small Automated CC Analyser	8
21 01 42		Medium Automated CC Analyser	9
21 01 43		Large Automated CC Analyser	
21 01 44		Extra large Automated CC Analyser	
21 01 51		Electrolyte - Flame photometer	10
21 01 52		Electrolyte/ISE only	
21 01 53		Electrolyte/ISE and other parameter	
21 02		Blood Gas Analyser	
21 02 01		Blood Gas (pO2 / pCO2 + pH)	
21 02 02		Blood Gas and Electrolytes	11
21 02 03		Blood Gas and CO-Oximeter	
21 02 04		Blood Gas, Electrolytes and CO-Oximeter	
21 02 06		Blood Gas, Electrolytes, CO-Oximeter and Other Parameters	
21 02 07		CO-Oximeter	
21 02 09		Other Blood Gas Instruments	
21 03		Electrophoresis Equipment	12
21 03 01		Manual Electrophoresis	
21 03 02		Densitometers	
21 03 03		Automated integrated Electrophoresis Systems	
21 04		Chromatography Equipment	
21 05		Urine Analyser	13
21 05 01		Manual Urine Analyzer	
21 05 02		Small Automated Urine Analyzer	
21 05 03		Medium Automated Urine Analyzer	
21 05 04		Large Automated Urine Analyzer	
21 06		Self Testing Equipment (Clin.Chem.)	
21 06 01		Blood Glucose Meter	
21 06 02		Cholesterol Meter	
21 06 03		Lactate Meter	
21 06 11		Multiparameter CC Meter	
21 09		Other Clinical Chemistry Instruments	

EDMA CODE	Company	Instrument-Name	Page
22	IMMUNOCHEMISTRY INSTRUMENTS		14
22 01	RIA Equipment		
22 01 01	RIA Counter		
22 02	Protein Analyser		
22 02 01	Small Protein Analyser		
22 02 02	Large Protein Analyser		
22 03	Immunoassay Analyser		
22 03 01	Manual I.A. Instruments / Readers		
22 03 02	Semi-automated I.A. Systems		15
22 03 03	Automated I.A. Systems		16
22 03 04	Batch / High throughput I.A. Systems		
22 03 05	Nucleic Acid Processors		
22 09	Other Immunochemistry Instruments		
23	HAEMATOLOGY INSTRUMENTS		17
23 01	Haemoglobin Analyser		
23 01 01	Dedicated Total Haemoglobin Analyser		
23 02	Blood Cell Counter		
23 02 01	No differential Blood Cell Counter		
23 02 02	2 / 3 part differential Blood Cell Counter		18
23 02 03	5 part differential low / medium end Blood Cell Counter		
23 02 04	5 part differential high end Blood Cell Counter		19
23 02 11	Dedicated Reticulocyte Analyser		
23 03	Flow Cytometer		
23 03 01	Flow Cytometer without Cell sorter		
23 03 02	Flow Cytometer with Cell sorter		
23 04	Coagulation Systems		20
23 04 01	Manual Coagulation Instruments		
23 04 02	Semi-automated Coagulation Systems		
23 04 03	Fully-automated Coagulation Systems		21
23 04 04	Aggregometer		
23 04 05	Selftesting Equipment (Coagulation)		
23 05	Blood Grouping Instruments		
23 05 01	Blood Grouping Instruments		
23 06	Histology Equipment (excl.microscopes)		
23 06 01	General Histology Instruments		
23 06 02	Microtomes		
23 09	Other Haematology Instruments		
23 09 01	Sediment rate Instruments		

EDMA CODE	Company	Instrument-Name	Page
24		MICROBIOLOGY INSTRUMENTS	22
24 01		Bacteriology Instruments	
24 01 01		Plate pouring Instruments	
24 01 02		Susceptibility and / or ID-systems	
24 01 03		Blood Culture Systems	
24 01 04		Dedicated Urine Screening MB Instruments	23
24 09		Other Microbiology Instruments	
26		OTHER CLINICAL INSTRUMENTS	
26 01		Clinical Centrifuges	
26 02		Sample Processors	
26 02 01		Clinical Chemistry Sample Processors	
26 02 02		Immuno Assays Sample Processors	
26 02 03		Haematology Sample Processors	
26 02 04		Microbiology Sample Processors	
26 02 09		Multiple Purpose Sample Processors	
26 03		Washers	
26 03 02		Immuno Assays Washers	
26 03 03		Haematology Washers	
26 03 09		Other Washers	
26 09		Other Other Clinical Instruments	
27		DATA MANAGEMENT SYSTEMS	
27 01		Data Management Hardware	
27 02		Data Management Software/Consumables	

EDMA CODE	Company	Instrument-Name
<u>Instrument Grouping</u>		
2	INSTRUMENTS	
21	CLINICAL CHEMISTRY INSTRUMENTS	
21 01	Chemistry Analyser	
21 01 01	Manual Photometer and Spectrophotometer	
	Bayer	ML 1 ML2 ML3
	Beckman	Model 42 DU 5
	BioMérieux	UV-VIS Jr UV-VIS 1 UV-VIS 2 VISUAL
	BioSystems	BTS 302 BTS 305 BTS 310
	Boehringer Mannheim	Photometer 4010 Photometer 4020 Photometer 5010
	Eppendorf	PCP 6121 ECOM 6122 ECOM 6124 F
	Gilford	103 203
	Hycel	ISIS
	Lange (Dr. Lange)	LP 1 LP 3 LP 6A LP 400 LP 500 Miniphotometer 8 Miniphotometer 20
	Merck	Eclipse Eclipse Plus
	Milton-Roy	Spectronic 1001+ Spectronic 1101 Spectronic 1201 Spectronic 20 / 20 D Spectronic 21 DV Spectronic 21 MV Spectronic 21 UVD Spectronic 3000 Spectronic 301 Spectronic 401 Spectronic 501 Spectronic 601
	Perkin-Elmer	Lambda 11 Lambda 12 Lambda 14
	Riele	PMC 2000 PM 210 PM 310
	Unicam	PU 8610 PU 8620
	Vitatron	IFP ISP

EDMA CODE	Company	Instrument-Name
21 01 11	Dedicated Glucose Systems	
	Beckman	Glucose-Analyzer 2
	Eppendorf	EBIO Compact
	Radiometer	ABL System 625
	Yellow Springs	YSI
21 01 12	Dedicated Glycated Haemoglobin /HbA1 Systems	
	ARKRAY KDK	Hi-Auto A1c HA 8140
	Bayer	DCA 2000
	Bio-Rad	DIAMAT
		DiaSTAT
		VARIANT
	Chiron	Glycomat
	Kontron	NN
	Menarini	HA 8110
		HA 8121
	HA 8140	
	Merck	Hitachi L 9000
21 01 19	Other Dedicated Single-Parameter CC Systems	
	Analox Instruments	LM10
	Beckman	Creatinine-Analyzer 2
		Urea-Analyzer 2
	Elvilogos	Bilirubinometer
21 01 29	Dedicated Multi-Parameter CC Systems	
	Analox Instruments	P-GM7 micro-Stat
	Eppendorf	EBIO Plus
21 01 31	Small Manual and Semi-automated CC Systems	
	Abbott	Vision
	Bayer	Seralyzer
	Boehringer Mannheim	Reflotron
		Reflotron Sprint
	Cholestech	L.D.X
	Ortho Diag. Systems (J&J)	Vitros DT 60/Vitros DT60
	Eppendorf	SPOTCHEM
	Menarini	SPOTCHEM
	Merck	EASY ST
	Roche	COBAS Ready (=SPOTCHEM)
21 01 32	Semi-automated CC Analyser	
	Beckman	Clinical System 700
	BPC BioSED	NN
	Eppendorf	EPAC 6140

EDMA CODE	Company	Instrument-Name
21 01 41	Small Automated CC Analyser	
	Abbott	ABA 100 / 200 VP Spectrum
	Bayer	RA 50 RA 500 RA 1000 RA XT
	Beckman	AS 4 / AS 8 Synchron CX3
	Boehringer Mannheim	Autolab Hitachi 704 Hitachi 705 Hitachi 902 Hitachi 904
	Chimica Diagnostica	Naxos
	Chiron	EXPRESS
	Dade (Dupont)	Dimension 380 Dimension ES
	Diagnostic Testing	Piccolo
	Eppendorf	ACP 5040 EASY EPOS ELAN
	Hycel	LISA 200 LISA 300
	Instrumentation Laboratory	Monarch 1000
	Kone	Progress Progress Plus Specific Specific Supra Pro
	Menarini	Classic/Artax
	Merck	Selectra
	Olympus	AU 200
	Roche	COBAS BIO COBAS MIRA COBAS MIRA S COBAS MIRA PLUS COBAS FARA COBAS FARA II
	Schiaparelli	ACE
	Sigma	SD2000
	Wako	20R

EDMA CODE	Company	Instrument-Name	
21 01 42	Medium Automated CC Analyser		
	Abbott	EPx	
	Bayer	RA2000 OPERA	
	Beckman	Synchron CX 3 DELTA Synchron CX 4 Synchron CX 4 DELTA Synchron CX 5 Synchron CX 5 DELTA	
	Boehringer Mannheim	Hitachi 911	
	Coulter	Optichem 500	
	Eppendorf	ERIS ETEC	
	Instrumentation Laboratory	Monarch PLUS Monarch 2000	
	Kone	Delta Dynamic Optima Ultra	
	Hycel	LISA 400	
	Menarini	Falcor 300 Falcor 600	
	Nova	Nucleus	
	Olympus	AU 510 AU 560	
	Ortho Diag. Systems (J&J)	Vitros 250 Vitros 500	
	Wako	30R	
	21 01 43	Large Automated CC Analyser	
		American Monitor	Perspective/Excel
		Bayer	AXON DAX 24 DAX 48
		Beckman	Synchron CX 7 / CX7 Delta
		Boehringer Mannheim	Hitachi 717 Hitachi 917
Dade (Dupont)		Paramax-RX Dimension 760 Dimension AR Dimension XL	
Instrumentation Laboratory		ILab 900	
Merck		MEGA	
Olympus		AU 600	
Ortho Diag. Systems (J&J)		Vitros 750 XRC Vitros 950 IRC/Vitros 950	
Roche		COBAS INTEGRA	
21 01 44		Extra large Automated CC Analyser	
		Bayer	CHEM 1 DAX 72 DAX 96
		Boehringer Mannheim	HIT 737 HIT 747 Series
	Instrumentation Laboratory	ILab 1800	
	Olympus	AU 800 AU 1000 AU 5000 Series AU 5200 Series	

EDMA CODE	Company	Instrument-Name	
21 01 51	Electrolyte - Flame photometer		
	Hycel	PHF 104 PHF 106	
	Instrumentation Laboratory	IL 943	
21 01 52	Electrolyte/ISE only		
	AVL	9180 Electrolyte Analyser	
	Beckman	E2A	
	Hycel	EASYLYTE	
	Instrumentation Laboratory	IL 614, 644, 654	
YSI	EasyLyte		
21 01 53	Electrolyte/ISE and other parameter		
	Beckman	ELISE	
21 02	Blood Gas Analyser		
21 02 01	Blood Gas (pO₂ / pCO₂ + pH)		
	AVL	AVL 945	
		AVL 990	
		AVL 995	
		OPTI 1	
		Compact 1	
		Compact 2	
	Chiron	CCD 238	
		CCD 248	
		CCD 278	
		CCD 280	
		CCD 840	
		Hewlett Packard	i-STAT
	Instrumentation Laboratory	IL 1304	
		IL 1306	
		IL 1312	
		IL 1420	
		IL 1610	
		IL 1620	
		IL BG3	
		Nova	Stat Profile 3
		Radiometer	ABL 5
	ABL 30		
	ABL 300		
	ABL 330		
	ABL 500		
	ABL 510		
ABL 520			

EDMA CODE	Company	Instrument-Name	
21 02 02	Blood Gas and Electrolytes		
	Chiron	CCD 348 CCD 850 CCD 860	
	Diametrics	IRMA	
	Instrumentation Laboratory	IL 1400	
		IL 1405	
		IL 1430	
		IL 1630	
		IL 1640	
		IL 1650	
		GEM Premier	
		GEM STAT	
	GEM 6+		
	GEM Premier Plus		
	Nova Biomedical	Stat Profile 1,2,4,5,7	
		Stat Profile Plus 9,10	
Stat Profile Ultra			
Radiometer	ABL 4		
	ABL 505		
	ABL 620		
AVL	OMNI		
21 02 03	Blood Gas and CO-Oximeter		
21 02 04	Blood Gas, Electrolytes and CO-Oximeter		
	AVL	OMNI	
	Chiron	CCD 845 CCD 855	
	Instrumentation Laboratory	Synthesis 25 / 35	
Radiometer	ABL 620		
21 02 06	Blood Gas, Electrolytes, CO-Oximeter and Other Parameters		
	Chiron	CCD 865	
Nova	Stat Profile Ultra		
21 02 07	CO-Oximeter		
	AVL	AVL 912	
	Chiron	CCD 270	
		CCD 815	
	Instrumentation Laboratory	IL 482	
IL 682			
Radiometer	OSM 3		
21 02 08	Osmometer		
	Gonotec	Osomat 030/030D Osomat 050	
	Slamed	Osmometer 800 cl	
	Roebing	Micro-Osmometer	
	Vitech Scientific	Model 3900	
	Wescor	Vapro	
	21 02 09	Sweat Analyzer	
	Wescor	Sweat-Chek Analyzer	
	21 02 90	Other Blood Gas Instruments	

EDMA CODE	Company	Instrument-Name
21 03	Electrophoresis Equipment	
21 03 01	Manual Electrophoresis	
	Beckman	PARAGON
	Cosmo	FED 3200
	Helena	NN
	Sebia	NN
21 03 02	Densitometers	
	Beckman	APPRAISE
	Helena	NN
	Olympus	Fractoscan Fractoscan Junior
	Sebia	HYDRASYS
21 03 03	Automated integrated Electrophoresis Systems	
	Beckman	APPRAISE PARAGON CZE 2000
	Eppendorf (JOOKOO)	CTE 1000 CTE 5000
	Helena	Concept 2000 REP Cardio REP PRINCE-CE AUTOIFE
	Menarini	HS 7
	Olympus	HITE SYSTEM 310 HITE SYSTEM 320 HITE SYSTEM 620
	Sebia	HYDRASYS
21 04	Chromatography Equipment	

EDMA CODE	Company	Instrument-Name
21 05	Urine Analyser	
21 05 01	Manual Urine Analyser	
	Menarini	MiniAution
21 05 02	Small Automated Urine Analyser	
	Bayer	Clinitek 50 Clinitek 100
	Boehringer Mannheim	Urilux Minitron
	IMACO	URIMAC
21 05 03	Medium Automated Urine Analyser	
	Bayer	Clinitek 200 Clinitek 200 +
	Boehringer Mannheim	Miditron Miditron Junior Urotron RL 9
	Behring	RAPIMAT
21 05 04	Large Automated Urine Analyser	
	Bayer	Clinitek Atlas
	Boehringer Mannheim	Supertron
	IRIS	900UDx
	Menarini	SuperAution
	Sysmex	UF-100
	Wako	NN
	Yeong Deong	NN
21 06	Self Testing Equipment (Clin.Chem.)	
21 06 01	Blood Glucose Meter	
21 06 02	Cholesterol Meter	
21 06 03	Lactate Meter	
21 06 11	Multi-Parameter CC Meter	
21 09	Other Clinical Chemistry Instruments	

EDMA CODE	Company	Instrument-Name
22	IMMUNOCHEMISTRY INSTRUMENTS	
22 01	RIA Equipment	
22 01 01	RIA Counter	
22 02	Protein Analyser	
22 02 01	Small Protein Analyser	
	Beckman	ICS
	Behring	TurbiTimer
	Bio-Whittaker	Fiax
	Nycomed	Card Reader
	Orion Diagnostica	Turbox
		QuikRead 3
22 02 02	Large Protein Analyser	
	Bayer	DPA-1
	Behring	BN 100
		BNA
		BN II
	Beckman	Array 2.0
		Array 360
		IMMAGE
	Boehringer Mannheim	KeySys
	Sanofi Diagnostics Pasteur	QM 300
22 03	Immunoassay Analyser	
22 03 01	Manual I.A. Instruments / Readers	
	Abbott	Matrix
		Quantum
		Quantum II
	Behring	Elisa-Photometer
		Elisa-Reader 301
		Elisa-Reader 311
		Berilux 250
		Microtrack Reader SE System
		ECOM 6123 D
	Biochem Immunosystems	Serazyme I
		Serazyme II
	BioMerieux	Axia Microreader
		Axia 2 Reader
	Boehringer Mannheim	ES 11
		ES 22
		ES 33
	Chiron	Magic Lite I
	Dynatech	System 5000
	Eppendorf	ECOM 6123 D
	Hybritech	Photon II
	Merck	Eclair
		LM-8
	Murex	Micro ELISA Reader 7000
		Micro ELISA MRX Reader
	Organon Teknika	Microelisa Reader 100
		Microelisa System
	Ortho Diagnostic Systems	Auto Reader 3
	Sanofi Diagnostics Pasteur	LP 400

EDMA CODE	Company	Instrument-Name
22 03 02	Semi-automated I.A. Systems	
	Abbott	ADx IMx / IMx Select TDx / TDxFLx
	Behring	BeriLux 400 BEP II OPUS OPUS PLUS ETS Plus Solaris
	Biochem Immunosystems	SR 1
	BioMerieux	VIDAS Mini-VIDAS
	Boehringer Mannheim	ES 300 ElecSys 1010
	Byk-Sangtek	LIA-mat S 300
	Chiron	Magic Lite II Magic Lite II+ Magic Lite SQ Allergy System
	Dade	Stratus I Stratus II
	DPC	Milenia Max Immulite
	Dade (Dupont)	ACA III ACA V ACA SX ACA STAR
	Hybritech	Photon Era
	Ortho Diag. Systems (J&J)	Amerlite
	Pharmacia	UniCAP
	Radim	BRIO
	Roche	COBAS EIA System
	Sorin	ETI-SYSTEM
	Tosoh	AIA 600
	Wallac / Pharmacia	Delfia

EDMA CODE	Company	Instrument-Name	
22 03 03	Automated I.A. Systems		
	Abbott	AxSYM	
	Bayer	Immuno-I	
	Behring	OPUS Magnum	
	Biochem Immunosystems	Labotec	
	Boehringer Mannheim	ES 600	
		ES 700	
		ElecSys 2010	
	Chiron	ACS 180	
		ACS 180 PLUS	
		Centaur	
	CIS-Bio	Cispack 4200	
	DPC	Immulite 2000	
	Hycor	HY-TEC 288	
	Kone	Quartus	
	Lisabio	Ligand Immuno System Assay	
	Merck (Biotrol)	MAGIA 7000	
		MAGIA 8000	
	Nichols Institute Diagnostic	Advantage	
	Organon Teknika	Tektime	
	Ortho Diag. Systems (J&J)	Amerlite PC	
		Vitros ECI	
	Radim	Ario	
	Roche	COBAS CORE	
		COBAS CORE II	
	Rosys	Plato 1300	
	Sanofi Diagnostics Pasteur	Access	
	Sorin	ETI-LAB	
	Tosoh	AIA 1200	
	Wallac / Pharmacia	AutoDELFIA	
		Unicap 100	
	22 03 04	Batch / High throughput I.A. Systems	
		Abbott	Commander / FPC PRISM
Beckman		Image	
Behring		BEP III	
		BEP XL (Micro Trak XL)	
Bio-Rad		Radius	
BIO-TEK		ELx808	
DPC / Biermann		MAX	
Dynatech		System 7000	
Hamilton		Microlab F.A.M.E.	
Hycor		HY TEC EIA System	
Merck		MIOS	
Murex		DIAS System	
		MIAS System	
Organon Teknika		MPP 3000	
Ortho Diagnostic Systems		Summit	
Rosys		Plato 3300	
Sanofi Diagnostics Pasteur		AMS/LP 400	
		Sonata	
Sorin		Omni System	
22 03 05		Nucleic Acid Processors	
		Abbott	LCx
		Becton Dickinson	ProbeTec
	Organon Teknika	NASBA QR System	
	Perkin Elmer	GenAmp PCR system 2400	
Roche	COBAS Amplicor		
22 09	Other Immunochemistry Instruments		

EDMA CODE	Company	Instrument-Name
23	HAEMATOLOGY INSTRUMENTS	
23 01	Haemoglobin Analyser	
23 01 01	Dedicated Total Haemoglobin Analyser	
23 02	Blood Cell Counter	
23 02 01	No differential Blood Cell Counter	
	Abbott	CELL-DYN 300 CELL-DYN 400 CELL-DYN 500 CELL-DYN 610 CELL-DYN 700 CELL-DYN 800 CELL-DYN 900
	Becton-Dickinson	QBC I QBC II QBC Autoread
	Biochemical Systems	Hema-Check 5
	Biochem Immunosystems	130 150 170 7500 8000
	Coulter	CBC-5 M430 MD 8 ZM M530
	Kone	AC 900 AC 920
	Nihon Kohden	Celltak MEK-5100
	Sysmex	F-520 F-820

EDMA CODE	Company	Instrument-Name	
23 02 02	2 / 3 part differential Blood Cell Counter		
	Abbott	CELL-DYN 1300 CELL-DYN 1400 CELL-DYN 1600 CELL-DYN 1700	
	Bayer	H1JR	
	Coulter	T540 T660 MD 16 MD II S770	
	Diatron	Twincell	
	Hycel	CELLY HEMACELL	
	Kone	AC 970	
	Nihon Kohden	Celltak Mini Celltak-Auto	
	Roche	COBAS MICROS COBAS MINOS ST COBAS MINOS STE COBAS MINOS STX COBAS MINOS STEX COBAS MINOS VET	
	Biochem Immunosystems	9000 9100 + 9118 + 9120 +	
	SEAC	GENIUS	
	Sysmex	K-800 K-1000	
	23 02 03	5 part differential low / medium end Blood Cell Counter	
		Abbott	CELL-DYN 2000
		Bayer	H1
		Biochemical Systems	Hema-Check 13
		Coulter	JS JR T890 S880 Onyx
		Hycel	HEMACELL PLUS
		Roche	COBAS HELIOS LMG COBAS HELIOS 5 Diff.
		Sysmex	NE-1500 K-4500

EDMA CODE	Company	Instrument-Name	
23 02 04	5 part differential high end Blood Cell Counter		
	Abbott	CELL-DYN 3000 CELL-DYN 3500 CELL-DYN 4000	
	Bayer	H 1 E H 2 H 3 RTC H 3 RTX	
	Coulter	STKR STKS MAXM ST	
	Lisabio	Gen S	
	Roche	Hemacell/Hemacell Plus COBAS ARGOS LMG COBAS ARGOS 5 Diff. COBAS VEGA	
	Swelab	AC920	
	Sysmex	SF-3000 SE-9000 NE-8000	
	23 02 11	Dedicated Reticulocyte Analyser	
		Sysmex	R-3000
		Sediment rate analyser	
		COFIM	SR 50
	Eurotec Medical System	VT-50/VT-75/VT-100	
	ALIFAX	Test 1	
23 03	Flow Cytometer		
23 03 01	Flow Cytometer without Cell sorter		
	Becton-Dickinson	FACScan FACStrak FACSCalibur FACSCount	
	Bio-Rad	NN	
	Coulter	Epics Profile II	
	Ortho Diag. Systems (J&J)	Cytonon Absolute	
	Partec	CA-II	
	Sysmex	R-2000/3000	
	23 03 02	Flow Cytometer with Cell sorter	
Becton-Dickinson		FACStarPlus FACSVantage FACSort	
Coulter		Epics Elite	
Partec	PAS-II		

EDMA CODE	Company	Instrument-Name
23 04	Coagulation Systems	
23 04 01	Manual Coagulation Instruments	
	Amelung	KC 4 KC 10 KC 10A
	Behnk	CL4 CL8
	Behring	Fibrintimer II Fibrintimer 10 Chromotimer
	BioMerieux	Option 2-8
	Boehringer Mannheim	Coaguheck/Coaguheck Plus
	ERMA	SSC-100
	Helena	Cascade M-4
	Labor	COA Data 1000-3000 COA System COA Screener
	Labsystems	FP 900 FP 910
	LRE	Coatron F Coatron F II
	Nycomed	Thrombotrack
	Organon Teknika	Coag-a-Mate XM
	Ortho Diagnostic Systems (MLA)	Electra 700 Koagulab CTS
	Roche	Fibrometer
	Stago	STA 2 STA 4 STA 8 STA 888
	Teco Medical Instruments	Coachrome IV
23 04 02	Semi-automated Coagulation Systems	
	Biodata	MCA 110
	Behnk	Thrombolyzer Chrom
	Boehringer Mannheim / Stago	Coasys 110
	Helena	Cascade 480
	Instrumentation Laboratory	ACL 100 ACL 200 ACL 300
	ITC-International Technidyne	ProTime Microcoagulation System
	LRE	Coatron Jr. Coatron Fx
	Organon Teknika	Coag-A-Mate X2 Coag-A-Mate RA 4
	Ortho Diagnostic Systems (MLA)	Electra 800 Electra 900 Electra 900 C Koagulab 16-605

EDMA CODE	Company	Instrument-Name
23 04 03	Fully-automated Coagulation Systems	
	Amelung	KC 40
	Behnk	Thrombolyser Thromboyzer Combi
	Behring	BFA BCT BCS
	BioMerieux	Hemolab
	Boehringer Mannheim / Stago	STA STA Compact
	Chiron	Biotrack 512
	Instrumentation Laboratory	ACL 1000 ACL 2000 ACL 3000 ACL 3000p ACL 6000 ACL Futura
	Organon Teknika	Coag-A-Mate MDA Coag-A-Mate MTX
	Ortho Diagnostic Systems (MLA)	Electra 1000 C Electra 1400 C Electra 1600 C
	Stratec	Stracomat
	Sysmex / Dade	CA-1000 CA-5000 CA-6000
23 04 04	Aggregometer	
	Biodata	PAP 4
	Ortho Diag. Systems / Dade	PFA-100
23 04 05	Selftesting Equipment (Coagulation)	
23 05	Blood Grouping Instruments	
23 05 01	Blood Grouping Instruments	
	Lifescreeen	QUASAR 240 QUASAR 400
	Menarini	BOSS
	Olympus	PK 80 PK 7100 PK 7200
	Ortho Diag. Systems (J&J)	AutoVue
23 06	Histology Equipment (excl.microscopes)	
23 06 01	General Histology Instruments	
	Pap Smear Screening NEOPATH	AutoPap 300 QC System
23 06 02	Microtomes	
23 09	Other Haematology Instruments	
23 09 01	Sediment rate Instruments	

EDMA CODE	Company	Instrument-Name
24	MICROBIOLOGY INSTRUMENTS	
24 01	Bacteriology Instruments	
24 01 01	Plate pouring Instruments	
24 01 02	Susceptibility and / or ID-systems	
	Dade	Touchscan Autoscan 4 Microscan Walkaway 40 Microscan Walkaway 96
	Becton Dickinson	Sceptor Autosceptor Sceptorpette
	BioMerieux	Vitek Junior Vitek 32 Vitek 60 Vitek 120 Vitek 240 mini API Uniscept Aladin ATB Expression
	Biokit	Videobact
	Difco	Pasco
	I2R	STRSCAN
	Mast	Radius Mastscan
	Radiometer	Sensititre
24 01 03	Blood Culture Systems	
	Becton Dickinson	Bactec 460 Tb Bactec 660 Bactec 730 Bactec 860 Bactec 9000 MB Bactec 9050 Bactec 9120 Bactec 9240
	BioMerieux	Vital 200 Vital 1200
	Difco	Sentinel ESP
	Organon Teknika	BacT/Alert 120 BacT/Alert 240
	Sanofi Diagnostics Pasteur	Bio Argos

EDMA CODE	Company	Instrument-Name
24 01 04	Dedicated Urine Screening MB Instruments	
	BioMerieux	Bact-T-Screen 2000 URIMAT
	Coral Biotechnology	Coral Wave 120
	Difco	Questor
	Lumac	Biocounter M2300 Biocounter M3000
	Orbec	Ramus
	Sysmex	UF-100
24 09	Other Microbiology Instruments	
26	OTHER CLINICAL INSTRUMENTS	
26 01	Clinical Centrifuges	
26 02	Sample Processors	
26 02 01	Clinical Chemistry Sample Processors	
26 02 02	Immuno Assays Sample Processors	
	Abbott	PPC
	Beckman	BIOMEK 2000
	Ortho Diag. Systems (J&J)	Summit Plus 2 Summit Processor
	Sanofi Diagnostics Pasteur	AMSS
	Tecan	GENESIS RSP 200
26 02 03	Haematology Sample Processors	
26 02 04	Microbiology Sample Processors	
26 02 09	Multiple Purpose Sample Processors	
	Boehringer Mannheim	VS 250 PPV 1200
26 03	Washers	
26 03 02	Immuno Assays Washers	
26 03 03	Haematology Washers	
26 03 09	Other Washers	
26 09	Other Other Clinical Instruments	
27	DATA MANAGEMENT SYSTEMS	
27 01	Data Management Hardware	
27 02	Data Management Software/Consumables	

EDMA CODE

Company

Instrument-Name

Alphabetical Order by Company

21	01	31	Abbott	Vision
21	01	41	Abbott	ABA 100 / 200
21	01	41	Abbott	VP
21	01	41	Abbott	Spectrum
21	01	42	Abbott	EPx
22	03	01	Abbott	Matrix
22	03	01	Abbott	Quantum
22	03	01	Abbott	Quantum II
22	03	02	Abbott	ADx
22	03	02	Abbott	IMx / IMx Select
22	03	02	Abbott	TDx / TDxFLx
22	03	03	Abbott	AxSYM
22	03	04	Abbott	Commander / FPC
22	03	04	Abbott	PRISM
22	03	05	Abbott	LCx
23	02	01	Abbott	CELL-DYN 300
23	02	01	Abbott	CELL-DYN 400
23	02	01	Abbott	CELL-DYN 500
23	02	01	Abbott	CELL-DYN 610
23	02	01	Abbott	CELL-DYN 700
23	02	01	Abbott	CELL-DYN 800
23	02	01	Abbott	CELL-DYN 900
23	02	02	Abbott	CELL-DYN 1300
23	02	02	Abbott	CELL-DYN 1400
23	02	02	Abbott	CELL-DYN 1600
23	02	02	Abbott	CELL-DYN 1700
23	02	03	Abbott	CELL-DYN 2000
23	02	04	Abbott	CELL-DYN 3000
23	02	04	Abbott	CELL-DYN 3500
23	02	04	Abbott	CELL-DYN 4000
26	02	02	Abbott	PPC
23	04	01	Amelung	KC 4
23	04	01	Amelung	KC 10
23	04	01	Amelung	KC 10A
23	04	03	Amelung	KC 40
21	01	43	American Monitor	Perspective/Excel
21	02	01	AVL	AVL 945
21	02	01	AVL	AVL 990
21	02	01	AVL	AVL 995
21	02	01	AVL	OPTI 1
21	02	01	AVL	Compact 1
21	02	01	AVL	Compact 2
21	02	04	AVL	OMNI
21	02	07	AVL	AVL 912
21	01	01	Bayer	ML 1
21	01	01	Bayer	ML2
21	01	01	Bayer	ML3
21	01	12	Bayer	DCA 2000
21	01	31	Bayer	Seralyzer
21	01	42	Bayer	RA2000
21	01	42	Bayer	OPERA
21	01	43	Bayer	AXON

EDMA CODE	Company	Instrument-Name
21 01 43	Bayer	DAX 24
21 01 43	Bayer	DAX 48
21 01 44	Bayer	CHEM 1
21 01 44	Bayer	DAX 72
21 01 44	Bayer	DAX 96
21 05 02	Bayer	Clinitek 50
21 05 02	Bayer	Clinitek 100
21 05 03	Bayer	Clinitek 200
21 05 03	Bayer	Clinitek 200 +
21 05 04	Bayer	Clinitek Atlas
22 02 02	Bayer	DPA-1
22 03 03	Bayer	Immuno-I
23 02 02	Bayer	H1JR
23 02 03	Bayer	H1
23 02 04	Bayer	H 1 E
23 02 04	Bayer	H 2
23 02 04	Bayer	H 3 RTC
23 02 04	Bayer	H 3 RTX
21 01 41	Bayer	RA 50
21 01 41	Bayer	RA 500
21 01 41	Bayer	RA 1000
21 01 41	Bayer	RA XT
21 01 01	Beckman	Model 42
21 01 01	Beckman	DU 5
21 01 11	Beckman	Glucose-Analyzer 2
21 01 19	Beckman	Creatinine-Analyzer 2
21 01 19	Beckman	Urea-Analyzer 2
21 01 32	Beckman	Clinical System 700
21 01 41	Beckman	AS 4 / AS 8
21 01 41	Beckman	Synchron CX3
21 01 42	Beckman	Synchron CX 3 DELTA
21 01 42	Beckman	Synchron CX 4
21 01 42	Beckman	Synchron CX 4 DELTA
21 01 42	Beckman	Synchron CX 5
21 01 42	Beckman	Synchron CX 5 DELTA
21 01 43	Beckman	Synchron CX 7
21 01 43	Beckman	Synchron CX 7 DELTA
21 01 52	Beckman	E2A
21 01 53	Beckman	ELISE
21 03 01	Beckman	PARAGON
21 03 02	Beckman	APPRAISE
21 03 03	Beckman	PARAGON CZE 2000
22 02 01	Beckman	ICS
22 02 02	Beckman	Array 2.0
22 02 02	Beckman	Array 360
22 02 02	Beckman	IMAGE
26 02 02	Beckman	BIOMEK 2000
22 03 05	Becton Dickinson	ProbeTec
24 01 02	Becton Dickinson	Sceptor
24 01 02	Becton Dickinson	Sceptorpette
24 01 03	Becton Dickinson	Bactec 460 Tb
24 01 03	Becton Dickinson	Bactec 660
24 01 03	Becton Dickinson	Bactec 730
24 01 03	Becton Dickinson	Bactec 860

EDMA CODE	Company	Instrument-Name
24 01 03	Becton Dickinson	Bactec 9000 MB
24 01 03	Becton Dickinson	Bactec 9050
24 01 03	Becton Dickinson	Bactec 9120
24 01 03	Becton Dickinson	Bactec 9240
23 02 01	Becton-Dickinson	QBC I
23 02 01	Becton-Dickinson	QBC II
23 02 01	Becton-Dickinson	QBC Autoread
23 03 01	Becton-Dickinson	FACScan
23 03 01	Becton-Dickinson	FACStrak
23 03 01	Becton-Dickinson	FACSCalibur
23 03 01	Becton-Dickinson	FACSCount
23 03 02	Becton-Dickinson	FACStar ^{Plus}
23 03 02	Becton-Dickinson	FACSVantage
23 03 02	Becton-Dickinson	FACSort
23 04 01	Behnk	CL4
23 04 01	Behnk	CL8
23 04 03	Behnk	Thrombolyser
21 05 03	Behring	RAPIMAT
22 02 01	Behring	TurbiTimer
22 02 02	Behring	BN 100
22 02 02	Behring	BNA
22 02 02	Behring	BN II
22 03 01	Behring	Elisa-Photometer
22 03 01	Behring	Elisa-Reader 301
22 03 01	Behring	Elisa-Reader 311
22 03 01	Behring	Berilux 250
22 03 01	Behring	Microtrack Reader SE System
22 03 01	Behring	ECOM 6123 D
22 03 02	Behring	Berilux 400
22 03 02	Behring	BEP II
22 03 02	Behring	OPUS
22 03 02	Behring	OPUS PLUS
22 03 02	Behring	ETS Plus
22 03 02	Behring	Solaris
22 03 03	Behring	OPUS Magnum
22 03 04	Behring	BEP III
22 03 04	Behring	BEP XL (Micro Trak XL)
23 04 01	Behring	Fibrintimer II
23 04 01	Behring	Fibrintimer 10
23 04 01	Behring	Chromotimer
23 04 03	Behring	BFA
23 04 03	Behring	BCT
23 04 03	Behring	BCS
21 01 12	Bio-Rad	DIAMAT
21 01 12	Bio-Rad	VARIANT
22 03 04	Bio-Rad	Radius
22 03 04	BIO-TEK	ELx808
22 02 01	Bio-Whittaker	Fiax
22 03 01	Biochem Immunosystems	Serazyme I
22 03 01	Biochem Immunosystems	Serazyme II
22 03 02	Biochem Immunosystems	SR 1
22 03 03	Biochem Immunosystems	Labotec
23 02 01	Biochem Immunosystems	130
23 02 01	Biochem Immunosystems	150

EDMA CODE	Company	Instrument-Name
23 02 01	Biochem Immunosystems	170
23 02 01	Biochem Immunosystems	7500
23 02 01	Biochem Immunosystems	8000
23 02 02	Biochem Immunosystems	9000
23 02 02	Biochem Immunosystems	9100 +
23 02 02	Biochem Immunosystems	9118 +
23 02 02	Biochem Immunosystems	9120 +
23 04 02	Biodata	MCA 110
23 04 04	Biodata	PAP 4
24 01 02	Biokit	Videobact
22 03 01	BioMerieux	Axia Microreader
22 03 01	BioMerieux	Axia 2 Reader
22 03 02	BioMerieux	VIDAS
22 03 02	BioMerieux	Mini-VIDAS
23 04 01	BioMerieux	Option 2-8
23 04 03	BioMerieux	Hemolab
24 01 02	BioMerieux	Vitek Junior
24 01 02	BioMerieux	Vitek 32
24 01 02	BioMerieux	Vitek 60
24 01 02	BioMerieux	Vitek 120
24 01 02	BioMerieux	Vitek 240
24 01 02	BioMerieux	Uniscept
24 01 02	BioMerieux	Aladin
24 01 02	BioMerieux	ATB Expression
24 01 03	BioMerieux	Vital 200
24 01 03	BioMerieux	Vital 1200
24 01 04	BioMerieux	Bact-T-Screen 2000
24 01 04	BioMerieux	URIMAT
21 01 01	BioSystems	BTS 302
21 01 01	BioSystems	BTS 305
21 01 01	BioSystems	BTS 310
21 01 01	Boehringer Mannheim	Photometer 4010
21 01 01	Boehringer Mannheim	Photometer 4020
21 01 31	Boehringer Mannheim	Reflotron
21 01 31	Boehringer Mannheim	Reflotron Sprint
21 01 41	Boehringer Mannheim	Hitachi 704
21 01 41	Boehringer Mannheim	Hitachi 705
21 01 41	Boehringer Mannheim	Hitachi 904
21 01 42	Boehringer Mannheim	Hitachi 911
21 01 43	Boehringer Mannheim	Hitachi 717
21 01 43	Boehringer Mannheim	Hitachi 917
21 01 44	Boehringer Mannheim	HIT 737
21 01 44	Boehringer Mannheim	HIT 747 Series
21 05 02	Boehringer Mannheim	Urilux
21 05 02	Boehringer Mannheim	Minitron
21 05 03	Boehringer Mannheim	Miditron
21 05 03	Boehringer Mannheim	Miditron Junior
21 05 03	Boehringer Mannheim	Urotron RL 9
21 05 04	Boehringer Mannheim	Supertron
22 02 02	Boehringer Mannheim	KeySys
22 03 01	Boehringer Mannheim	ES 11
22 03 01	Boehringer Mannheim	ES 22
22 03 01	Boehringer Mannheim	ES 33
22 03 02	Boehringer Mannheim	ES 300

EDMA CODE	Company	Instrument-Name
22 03 02	Boehringer Mannheim	ElecSys 1010
22 03 03	Boehringer Mannheim	ES 600
22 03 03	Boehringer Mannheim	ES 700
22 03 03	Boehringer Mannheim	ElecSys 2010
26 02 09	Boehringer Mannheim	PPV
23 04 02	Boehringer Mannheim / Stago	Coasys 110
23 04 03	Boehringer Mannheim / Stago	STA
23 04 03	Boehringer Mannheim / Stago	STA Compact
22 03 02	Byk-Sangtek	LIA-mat S 300
21 01 12	Chiron (Ciba Corning)	Glycomat
21 01 41	Chiron (Ciba Corning)	EXPRESS
21 02 01	Chiron (Ciba Corning)	CCD 238
21 02 01	Chiron (Ciba Corning)	CCD 248
21 02 01	Chiron (Ciba Corning)	CCD 278
21 02 01	Chiron (Ciba Corning)	CCD 280
21 02 01	Chiron (Ciba Corning)	CCD 840
21 02 02	Chiron (Ciba Corning)	CCD 348
21 02 02	Chiron (Ciba Corning)	CCD 850
21 02 02	Chiron (Ciba Corning)	CCD 860
21 02 04	Chiron (Ciba Corning)	CCD 845
21 02 04	Chiron (Ciba Corning)	CCD 855
21 02 04	Chiron (Ciba Corning)	CCD 865
21 02 06	Chiron (Ciba Corning)	CCD 865
21 02 07	Chiron (Ciba Corning)	CCD 270
21 02 07	Chiron (Ciba Corning)	CCD 815
22 03 01	Chiron (Ciba Corning)	Magic Lite I
22 03 02	Chiron (Ciba Corning)	Magic Lite II
22 03 02	Chiron (Ciba Corning)	Magic Lite II+
22 03 02	Chiron (Ciba Corning)	Magic Lite SQ Allergy System
22 03 03	Chiron (Ciba Corning)	ACS 180
22 03 03	Chiron (Ciba Corning)	ACS 180 PLUS
23 04 03	Chiron (Ciba Corning)	Biotrack 512
21 01 31	Cholestech	L.D.X
22 03 03	CIS-Bio	Cispack 4200
24 01 04	Coral Biotechnology	Coral Wave 120
21 03 01	COSMO	FED 3200
21 01 42	Coulter	Optichem 500
23 02 01	Coulter	CBC-5
23 02 01	Coulter	M430
23 02 01	Coulter	MD 8
23 02 01	Coulter	ZM
23 02 01	Coulter	M530
23 02 02	Coulter	T540
23 02 02	Coulter	T660
23 02 02	Coulter	MD 16
23 02 02	Coulter	MD II
23 02 02	Coulter	S770
23 02 03	Coulter	JS
23 02 03	Coulter	JR
23 02 03	Coulter	T890
23 02 03	Coulter	S880
23 02 03	Coulter	Onyx
23 02 04	Coulter	STKR
23 02 04	Coulter	STKS

EDMA CODE	Company	Instrument-Name
23 02 04	Coulter	MAXM
23 02 04	Coulter	ST
23 03 01	Coulter	Epics Profile II
23 03 02	Coulter	Epics Elite
22 03 02	Dade	Stratus I
22 03 02	Dade	Stratus II
24 01 02	Dade	Touchscan
24 01 02	Dade	Autoscan 4
24 01 02	Dade	Microscan Walkaway 40
24 01 02	Dade	Microscan Walkaway 96
24 01 02	Dade	Autosceptor
21 01 41	Dade (Dupont)	Dimension 380
21 01 41	Dade (Dupont)	Dimension ES
21 01 43	Dade (Dupont)	Paramax-RX
21 01 43	Dade (Dupont)	Dimension 760
21 01 43	Dade (Dupont)	Dimension AR
21 01 43	Dade (Dupont)	Dimension XL
22 03 02	Dade (Dupont)	ACA III
22 03 02	Dade (Dupont)	ACA V
22 03 02	Dade (Dupont)	ACA SX
22 03 02	Dade (Dupont)	ACA STAR
21 02 02	Diametrics	IRMA
24 01 02	Difco	Pasco
24 01 03	Difco	Sentinel
24 01 03	Difco	ESP
24 01 04	Difco	Questor
22 03 02	DPC	Milenia Max
22 03 02	DPC	Immulite
22 03 04	DPC / Biermann	MAX
22 03 01	Dynatech	System 5000
22 03 04	Dynatech	System 7000
21 01 01	Eppendorf	PCP 6121
21 01 01	Eppendorf	ECOM 6122
21 01 01	Eppendorf	ECOM 6124 F
21 01 11	Eppendorf	EBIO Compact
21 01 29	Eppendorf	EBIO Plus
21 01 31	Eppendorf	SPOTCHEM
21 01 32	Eppendorf	EPAC 6140
21 01 41	Eppendorf	ACP 5040
21 01 41	Eppendorf	EASY
21 01 41	Eppendorf	EPOS
21 01 41	Eppendorf	ELAN
21 01 42	Eppendorf	ERIS
21 01 42	Eppendorf	ETEC
22 03 01	Eppendorf	ECOM 6123 D
21 03 03	Eppendorf (JOKOO)	CTE 1000
21 03 03	Eppendorf (JOKOO)	CTE 5000
21 01 01	Gilford	103
21 01 01	Gilford	203
22 03 04	Hamilton	Microlab F.A.M.E.
21 03 01	Helena	NN
21 03 02	Helena	NN
21 03 03	Helena	Concept 2000
21 03 03	Helena	REP

EDMA CODE	Company	Instrument-Name
21 03 03	Helena	Cardio REP
21 03 03	Helena	PRINCE-CE
23 04 01	Helena	Cascade M-4
23 04 02	Helena	Cascade 480
22 03 01	Hybritech	Photon II
22 03 02	Hybritech	Photon Era
24 01 02	I2R	STRSCAN
21 05 02	IMACO	URIMAC
21 01 41	Instrument. Lab. (IL)	Monarch 1000
21 01 42	Instrument. Lab. (IL)	Monarch PLUS
21 01 42	Instrument. Lab. (IL)	Monarch 2000
21 01 43	Instrument. Lab. (IL)	ILab 900
21 01 44	Instrument. Lab. (IL)	ILab 1800
21 01 51	Instrument. Lab. (IL)	IL 943
21 01 52	Instrument. Lab. (IL)	IL 614, 644, 654
21 02 01	Instrument. Lab. (IL)	IL 1304
21 02 01	Instrument. Lab. (IL)	IL 1306
21 02 01	Instrument. Lab. (IL)	IL 1312
21 02 01	Instrument. Lab. (IL)	IL 1420
21 02 01	Instrument. Lab. (IL)	IL 1610
21 02 01	Instrument. Lab. (IL)	IL 1620
21 02 01	Instrument. Lab. (IL)	IL BG3
21 02 02	Instrument. Lab. (IL)	IL 1400
21 02 02	Instrument. Lab. (IL)	IL 1405
21 02 02	Instrument. Lab. (IL)	IL 1430
21 02 02	Instrument. Lab. (IL)	IL 1630
21 02 02	Instrument. Lab. (IL)	IL 1640
21 02 02	Instrument. Lab. (IL)	IL 1650
21 02 02	Instrument. Lab. (IL)	GEM PREMIER
21 02 02	Instrument. Lab. (IL)	GEM STAT
21 02 02	Instrument. Lab. (IL)	GEM 6+
21 02 02	Instrument. Lab. (IL)	GEM PREMIER PLUS
21 02 04	Instrument. Lab. (IL)	Synthesis 25 /35
21 02 07	Instrument. Lab. (IL)	IL 482
21 02 07	Instrument. Lab. (IL)	IL 682
23 04 02	Instrument. Lab. (IL)	ACL 100
23 04 02	Instrument. Lab. (IL)	ACL 200
23 04 02	Instrument. Lab. (IL)	ACL 300
23 04 03	Instrument. Lab. (IL)	ACL 1000
23 04 03	Instrument. Lab. (IL)	ACL 2000
23 04 03	Instrument. Lab. (IL)	ACL 3000
23 04 03	Instrument. Lab. (IL)	ACL 3000p
23 04 03	Instrument. Lab. (IL)	ACL 6000
23 04 03	Instrument. Lab. (IL)	ACL Futura
21 01 41	Kone	Progress
21 01 41	Kone	Progress Plus
21 01 41	Kone	Specific
21 01 41	Kone	Specific Supra
21 01 42	Kone	Delta
21 01 42	Kone	Dynamic
21 01 42	Kone	Optima
21 01 42	Kone	Ultra
23 02 01	Kone	AC 900
23 02 01	Kone	AC 920

EDMA CODE	Company	Instrument-Name
23 02 02	Kone	AC 970
22 03 03	Kone	Quartus
23 04 01	Labor	COA Data 1000-3000
23 04 01	Labor	COA System
23 04 01	Labor	COA Screener
23 04 01	Labsystems	FP 900
23 04 01	Labsystems	FP 910
21 01 01	Lange (Dr. Lange)	LP 1
21 01 01	Lange (Dr. Lange)	LP 3
21 01 01	Lange (Dr. Lange)	LP 6A
21 01 01	Lange (Dr. Lange)	LP 400
21 01 01	Lange (Dr. Lange)	LP 500
21 01 01	Lange (Dr. Lange)	Miniphotometer 8
23 05 01	Lifescreeen	QUASAR 240
23 05 01	Lifescreeen	QUASAR 400
23 04 01	LRE	Coatron F
23 04 01	LRE	Coatron F II
23 04 02	LRE	Coatron Jr.
23 04 02	LRE	Coatron Fx
24 01 04	Lumac	Biocounter M2300
24 01 04	Lumac	Biocounter M3000
24 01 02	Mast	Radius
24 01 02	Mast	Mastscan
21 01 12	Menarini	HA 8110
21 01 12	Menarini	HA 8121
21 01 12	Menarini	HA 8140
21 01 31	Menarini	SPOTCHEM
21 01 41	Menarini	Classic/Artax
21 01 42	Menarini	Falcor 300
21 01 42	Menarini	Falcor 600
21 03 03	Menarini	HS 7
21 05 01	Menarini	MiniAution
21 05 04	Menarini	SuperAution
23 05 01	Menarini	BOSS
21 01 01	Merck	Eclipse
21 01 01	Merck	Eclipse Plus
21 01 12	Merck	Hitachi L 9000
21 01 31	Merck	EASY ST
21 01 41	Merck	Selectra
21 01 43	Merck	MEGA
22 03 01	Merck	Eclair
22 03 01	Merck	LM-8
22 03 04	Merck	MIOS
22 03 03	Merck (Biotrol)	MAGIA 7000
22 03 03	Merck (Biotrol)	MAGIA 8000
21 01 01	Milton-Roy	Spectronic 1001+
21 01 01	Milton-Roy	Spectronic 1101
21 01 01	Milton-Roy	Spectronic 1201
21 01 01	Milton-Roy	Spectronic 20 / 20 D
21 01 01	Milton-Roy	Spectronic 21 DV
21 01 01	Milton-Roy	Spectronic 21 MV
21 01 01	Milton-Roy	Spectronic 21 UVD
21 01 01	Milton-Roy	Spectronic 3000
21 01 01	Milton-Roy	Spectronic 301

EDMA CODE	Company	Instrument-Name
21 01 01	Milton-Roy	Spectronic 401
21 01 01	Milton-Roy	Spectronic 501
21 01 01	Milton-Roy	Spectronic 601
22 03 01	Murex	Micro ELISA Reader 7000
22 03 01	Murex	Micro ELISA MRX Reader
22 03 04	Murex	DIAS System
22 03 04	Murex	MIAS System
23 02 01	Nihon Kohden	Celltak MEK-5100
23 02 02	Nihon Kohden	Celltak-Auto
21 01 42	Nova	Nucleus
21 02 01	Nova	Stat Profile 3
21 02 02	Nova Biomedical	Stat Profile1,2,4,5,7
21 02 02	Nova Biomedical	Stat Profile Plus 9,10
22 02 01	Nycomed	Card Reader
23 04 01	Nycomed	Thrombotrack
21 01 41	Olympus	AU 200
21 01 42	Olympus	AU 510
21 01 42	Olympus	AU 560
21 01 43	Olympus	AU 600
21 01 44	Olympus	AU 800
21 01 44	Olympus	AU 1000
21 01 44	Olympus	AU 5000 Series
21 01 44	Olympus	AU 5200 Series
21 03 02	Olympus	Fractoscan
21 03 02	Olympus	Fractoscan Junior
21 03 03	Olympus	HITE SYSTEM 310
21 03 03	Olympus	HITE SYSTEM 320
21 03 03	Olympus	HITE SYSTEM 620
23 05 01	Olympus	PK 80
23 05 01	Olympus	PK 7100
23 05 01	Olympus	PK 7200
24 01 04	Orbec	Ramus
23 04 01	Organon Teknika	Coag-a-Mate XM
22 03 01	Organon Teknika	Microelisa Reader 100
22 03 01	Organon Teknika	Microelisa System
22 03 02	Organon Teknika	Auraflex
22 03 03	Organon Teknika	Tektime
22 03 04	Organon Teknika	MPP 3000
22 03 05	Organon Teknika	NASBA QR System
23 04 02	Organon Teknika	Coag-A-Mate X2
23 04 02	Organon Teknika	Coag-A-Mate RA 4
23 04 03	Organon Teknika	Coag-A-Mate MDA
23 04 03	Organon Teknika	Coag-A-Mate MTX
24 01 03	Organon Teknika	BacT/Alert 120
24 01 03	Organon Teknika	BacT/Alert 240
22 02 01	Orion Diagnostica	Turbox
22 02 01	Orion Diagnostica	QuikRead 3
21 01 31	Ortho Diag.Syst. (J&J)	VITROS(Ektachem) DT 60
21 01 42	Ortho Diag.Syst. (J&J)	VITROS(Ektachem) 250
21 01 42	Ortho Diag.Syst. (J&J)	VITROS(Ektachem) 500
21 01 43	Ortho Diag.Syst. (J&J)	VITROS(Ektachem) 750 XRC
21 01 43	Ortho Diag.Syst. (J&J)	VITROS(Ektachem) 950 IRC
22 03 02	Ortho Diag.Syst. (J&J)	Amerlite
22 03 03	Ortho Diag.Syst. (J&J)	Amerlite PC

EDMA CODE	Company	Instrument-Name
22 03 03	Ortho Diag.Syst. (J&J)	VITROS ECI
23 03 01	Ortho Diag.Syst. (J&J)	Cytoron Absolute
26 02 02	Ortho Diag.Syst. (J&J)	Summit Plus 2
26 02 02	Ortho Diag.Syst. (J&J)	Summit Processor
22 03 01	Ortho Diagnostic Systems	Auto Reader 3
22 03 04	Ortho Diagnostic Systems	Summit
23 04 01	Ortho Diagnostic Systems (MLA)	Electra 700
23 04 01	Ortho Diagnostic Systems (MLA)	Koagulab CTS
23 04 02	Ortho Diagnostic Systems (MLA)	Electra 800
23 04 02	Ortho Diagnostic Systems (MLA)	Electra 900
23 04 02	Ortho Diagnostic Systems (MLA)	Electra 900 C
23 04 02	Ortho Diagnostic Systems (MLA)	Koagulab 16-605
23 04 03	Ortho Diagnostic Systems (MLA)	Electra 1000 C
23 04 03	Ortho Diagnostic Systems (MLA)	Electra 1400 C
23 04 03	Ortho Diagnostic Systems (MLA)	Electra 1600 C
23 04 04	Ortho Diagnostic Systems / Dade	PFA-100
23 03 01	Partec	CA-II
23 03 02	Partec	PAS-II
21 01 01	Perkin-Elmer	Lambda 11
21 01 01	Perkin-Elmer	Lambda 12
21 01 01	Perkin-Elmer	Lambda 14
22 03 02	Pharmacia	UniCAP
22 03 02	Radim	BRIO
21 02 01	Radiometer	ABL 5
21 02 01	Radiometer	ABL 30
21 02 01	Radiometer	ABL 300
21 02 01	Radiometer	ABL 330
21 02 01	Radiometer	ABL 500
21 02 01	Radiometer	ABL 510
21 02 01	Radiometer	ABL 520
21 02 02	Radiometer	ABL 4
21 02 02	Radiometer	ABL 505
21 02 02	Radiometer	ABL 620
21 02 04	Radiometer	ABL 620
21 02 07	Radiometer	OSM 3
24 01 02	Radiometer	Sensititre
21 01 01	Riele	PMC 2000
21 01 01	Riele	PM 210
21 01 01	Riele	PM 310
21 01 31	Roche	COBAS Ready (=SPOTCHEM)
21 01 41	Roche	COBAS BIO
21 01 41	Roche	COBAS MIRA
21 01 41	Roche	COBAS MIRA S
21 01 41	Roche	COBAS MIRA PLUS
21 01 41	Roche	COBAS FARA
21 01 41	Roche	COBAS FARA II
21 01 43	Roche	COBAS INTEGRA
22 03 02	Roche	COBAS EIA System
22 03 03	Roche	COBAS CORE
22 03 03	Roche	COBAS CORE II
22 03 05	Roche	COBAS Amplicor
23 02 02	Roche	COBAS MICROS
23 02 02	Roche	COBAS MINOS ST
23 02 02	Roche	COBAS MINOS STE

EDMA CODE	Company	Instrument-Name
23 02 02	Roche	COBAS MINOS STX
23 02 02	Roche	COBAS MINOS STEX
23 02 02	Roche	COBAS MINOS VET
23 02 03	Roche	COBAS HELIOS LMG
23 02 03	Roche	COBAS HELIOS 5 Diff.
23 02 04	Roche	COBAS ARGOS LMG
23 02 04	Roche	COBAS ARGOS 5 Diff.
23 02 04	Roche	COBAS VEGA
23 04 01	Roche	Fibrometer
22 03 03	Rosys	Plato 1300
22 03 04	Rosys	Plato 3300
22 02 02	Sanofi Diagnostics Pasteur	QM 300
22 03 01	Sanofi Diagnostics Pasteur	LP 400
22 03 03	Sanofi Diagnostics Pasteur	Access
22 03 04	Sanofi Diagnostics Pasteur	AMS/LP 400
22 03 04	Sanofi Diagnostics Pasteur	Sonata
24 01 03	Sanofi Diagnostics Pasteur	Bio Argos
26 02 02	Sanofi Diagnostics Pasteur	AMSS
23 02 02	SEAC	GENIUS
21 03 01	Sebia	NN
21 03 02	Sebia	NN
21 03 03	Sebia	HYDRASYS
22 03 03	Sorin	ETI System
22 03 04	Sorin	Omni System
23 04 01	Stago	STA 2
23 04 01	Stago	STA 4
23 04 01	Stago	STA 8
23 04 01	Stago	STA 888
23 04 03	Stratec	Stracomat
23 02 01	Sysmex	F-520
23 02 01	Sysmex	F-820
23 02 02	Sysmex	K-800
23 02 02	Sysmex	K-1000
23 02 03	Sysmex	NE-1500
23 02 03	Sysmex	K-4500
23 02 04	Sysmex	SF-3000
23 02 04	Sysmex	SE-9000
23 02 04	Sysmex	NE-8000
23 02 11	Sysmex	R-3000
23 03 01	Sysmex	R-2000/3000
24 01 04	Sysmex	UF-100
23 04 03	Sysmex / Dade	CA-1000
23 04 03	Sysmex / Dade	CA-5000
23 04 03	Sysmex / Dade	CA-6000
26 02 02	Tecan	GENESIS RSP 200
22 03 02	Tosoh	AIA 600
22 03 03	Tosoh	AIA 1200
21 01 01	Unicam	PU 8610
21 01 01	Unicam	PU 8620
21 01 01	Vitatron	IFP
21 01 01	Vitatron	ISP
21 01 41	Wako	20R
21 05 04	Wako	NN
21 01 42	Wako	30R

EDMA CODE	Company	Instrument-Name
22 03 02	Wallac / Pharmacia	Delfia
22 03 03	Wallac / Pharmacia	AutoDELFIA
21 01 11	Yellow Springs	YSI
21 05 04	Yeong Deong	NN

EDMA CODE

Company

Instrument-Name

Alphabetical Order by Instrument Name

21 01 01	Gilford	103
23 02 01	Biochem Immunosystems	130
23 02 01	Biochem Immunosystems	150
23 02 01	Biochem Immunosystems	170
21 01 01	Gilford	203
23 02 01	Biochem Immunosystems	7500
23 02 01	Biochem Immunosystems	8000
23 02 02	Biochem Immunosystems	9000
21 01 41	Wako	20R
21 01 42	Wako	30R
23 02 02	Biochem Immunosystems	9100 +
23 02 02	Biochem Immunosystems	9118 +
23 02 02	Biochem Immunosystems	9120 +
21 01 41	Abbott	ABA 100 / 200
21 02 01	Radiometer	ABL 30
21 02 01	Radiometer	ABL 300
21 02 01	Radiometer	ABL 330
21 02 02	Radiometer	ABL 4
21 02 01	Radiometer	ABL 5
21 02 01	Radiometer	ABL 500
21 02 02	Radiometer	ABL 505
21 02 01	Radiometer	ABL 510
21 02 01	Radiometer	ABL 520
21 02 02	Radiometer	ABL 620
21 02 04	Radiometer	ABL 620
23 02 01	Kone	AC 900
23 02 01	Kone	AC 920
23 02 02	Kone	AC 970
22 03 02	Dade (Dupont)	ACA III
22 03 02	Dade (Dupont)	ACA STAR
22 03 02	Dade (Dupont)	ACA SX
22 03 02	Dade (Dupont)	ACA V
22 03 03	Sanofi Diagnostics Pasteur	Access
23 04 02	Instrument. Lab. (IL)	ACL 100
23 04 03	Instrument. Lab. (IL)	ACL 1000
23 04 02	Instrument. Lab. (IL)	ACL 200
23 04 03	Instrument. Lab. (IL)	ACL 2000
23 04 02	Instrument. Lab. (IL)	ACL 300
23 04 03	Instrument. Lab. (IL)	ACL 3000
23 04 03	Instrument. Lab. (IL)	ACL 3000p
23 04 03	Instrument. Lab. (IL)	ACL 6000
23 04 03	Instrument. Lab. (IL)	ACL Futura
21 01 41	Eppendorf	ACP 5040
22 03 03	Chiron (Ciba Corning)	ACS 180
22 03 03	Chiron (Ciba Corning)	ACS 180 PLUS
22 03 02	Abbott	ADx
22 03 03	Tosoh	AIA 1200
22 03 02	Tosoh	AIA 600
24 01 02	BioMerieux	Aladin
22 03 02	Ortho Diag.Syst. (J&J)	Amerlite
22 03 03	Ortho Diag.Syst. (J&J)	Amerlite PC
22 03 04	Sanofi Diagnostics Pasteur	AMS/LP 400

EDMA CODE	Company	Instrument-Name
26 02 02	Sanofi Diagnostics Pasteur	AMSS
21 03 02	Beckman	APPRAISE
22 02 02	Beckman	Array 2.0
22 02 02	Beckman	Array 360
21 01 41	Beckman	AS 4 / AS 8
24 01 02	BioMerieux	ATB Expression
21 01 44	Olympus	AU 1000
21 01 41	Olympus	AU 200
21 01 44	Olympus	AU 5000 Series
21 01 42	Olympus	AU 510
21 01 44	Olympus	AU 5200 Series
21 01 42	Olympus	AU 560
21 01 43	Olympus	AU 600
21 01 44	Olympus	AU 800
22 03 02	Organon Teknika	Auraflex
22 03 01	Ortho Diagnostic Systems	Auto Reader 3
22 03 03	Wallac / Pharmacia	AutoDELFLIA
24 01 02	Dade	Autoscan 4
24 01 02	Dade	Autosceptor
21 02 07	AVL	AVL 912
21 02 01	AVL	AVL 945
21 02 01	AVL	AVL 990
21 02 01	AVL	AVL 995
22 03 01	BioMerieux	Axia 2 Reader
22 03 01	BioMerieux	Axia Microreader
21 01 43	Bayer	AXON
22 03 03	Abbott	AxSYM
24 01 04	BioMerieux	Bact-T-Screen 2000
24 01 03	Organon Teknika	BacT/Alert 120
24 01 03	Organon Teknika	BacT/Alert 240
24 01 03	Becton Dickinson	Bactec 460 Tb
24 01 03	Becton Dickinson	Bactec 660
24 01 03	Becton Dickinson	Bactec 730
24 01 03	Becton Dickinson	Bactec 860
24 01 03	Becton Dickinson	Bactec 9000 MB
24 01 03	Becton Dickinson	Bactec 9050
24 01 03	Becton Dickinson	Bactec 9120
24 01 03	Becton Dickinson	Bactec 9240
23 04 03	Behring	BCS
23 04 03	Behring	BCT
22 03 02	Behring	BEP II
22 03 04	Behring	BEP III
22 03 04	Behring	BEP XL (Micro Trak XL)
22 03 01	Behring	Berilux 250
22 03 02	Behring	Berilux 400
23 04 03	Behring	BFA
24 01 03	Sanofi Diagnostics Pasteur	Bio Argos
24 01 04	Lumac	Biocounter M2300
24 01 04	Lumac	Biocounter M3000
26 02 02	Beckman	BIOMEK 2000
23 04 03	Chiron (Ciba Corning)	Biotrack 512
22 02 02	Behring	BN 100
22 02 02	Behring	BN II
22 02 02	Behring	BNA

EDMA CODE	Company	Instrument-Name
23 05 01	Menarini	BOSS
22 03 02	Radim	BRIO
21 01 01	BioSystems	BTS 302
21 01 01	BioSystems	BTS 305
21 01 01	BioSystems	BTS 310
23 04 03	Sysmex / Dade	CA-1000
23 04 03	Sysmex / Dade	CA-5000
23 04 03	Sysmex / Dade	CA-6000
23 03 01	Partec	CA-II
22 02 01	Nycomed	Card Reader
21 03 03	Helena	Cardio REP
23 04 02	Helena	Cascade 480
23 04 01	Helena	Cascade M-4
23 02 01	Coulter	CBC-5
21 02 01	Chiron (Ciba Corning)	CCD 238
21 02 01	Chiron (Ciba Corning)	CCD 248
21 02 07	Chiron (Ciba Corning)	CCD 270
21 02 01	Chiron (Ciba Corning)	CCD 278
21 02 01	Chiron (Ciba Corning)	CCD 280
21 02 02	Chiron (Ciba Corning)	CCD 348
21 02 07	Chiron (Ciba Corning)	CCD 815
21 02 01	Chiron (Ciba Corning)	CCD 840
21 02 04	Chiron (Ciba Corning)	CCD 845
21 02 02	Chiron (Ciba Corning)	CCD 850
21 02 04	Chiron (Ciba Corning)	CCD 855
21 02 02	Chiron (Ciba Corning)	CCD 860
21 02 04	Chiron (Ciba Corning)	CCD 865
21 02 06	Chiron (Ciba Corning)	CCD 865
23 02 02	Abbott	CELL-DYN 1300
23 02 02	Abbott	CELL-DYN 1400
23 02 02	Abbott	CELL-DYN 1600
23 02 02	Abbott	CELL-DYN 1700
23 02 03	Abbott	CELL-DYN 2000
23 02 01	Abbott	CELL-DYN 300
23 02 04	Abbott	CELL-DYN 3000
23 02 04	Abbott	CELL-DYN 3500
23 02 01	Abbott	CELL-DYN 400
23 02 04	Abbott	CELL-DYN 4000
23 02 01	Abbott	CELL-DYN 500
23 02 01	Abbott	CELL-DYN 610
23 02 01	Abbott	CELL-DYN 700
23 02 01	Abbott	CELL-DYN 800
23 02 01	Abbott	CELL-DYN 900
23 02 01	Nihon Kohden	Celltak MEK-5100
23 02 02	Nihon Kohden	Celltak-Auto
21 01 44	Bayer	CHEM 1
23 04 01	Behring	Chromotimer
22 03 03	CIS-Bio	Cispack 4200
23 04 01	Behnk	CL4
23 04 01	Behnk	CL8
21 01 41	Menarini	Classic/Artax
21 01 32	Beckman	Clinical System 700
21 05 02	Bayer	Clinitek 100
21 05 03	Bayer	Clinitek 200

EDMA CODE	Company	Instrument-Name
21 05 03	Bayer	Clinitek 200 +
21 05 02	Bayer	Clinitek 50
21 05 04	Bayer	Clinitek Atlas
23 04 01	Labor	COA Data 1000-3000
23 04 01	Labor	COA Screener
23 04 01	Labor	COA System
23 04 03	Organon Teknika	Coag-A-Mate MDA
23 04 03	Organon Teknika	Coag-A-Mate MTX
23 04 02	Organon Teknika	Coag-A-Mate RA 4
23 04 02	Organon Teknika	Coag-A-Mate X2
23 04 01	Organon Teknika	Coag-a-Mate XM
23 04 02	Boehringer Mannheim / Stago	Coasys 110
23 04 01	LRE	Coatron F
23 04 01	LRE	Coatron F II
23 04 02	LRE	Coatron Fx
23 04 02	LRE	Coatron Jr.
22 03 05	Roche	COBAS Amplicor
23 02 04	Roche	COBAS ARGOS 5 Diff.
23 02 04	Roche	COBAS ARGOS LMG
21 01 41	Roche	COBAS BIO
22 03 03	Roche	COBAS CORE
22 03 03	Roche	COBAS CORE II
22 03 02	Roche	COBAS EIA System
21 01 41	Roche	COBAS FARA
21 01 41	Roche	COBAS FARA II
23 02 03	Roche	COBAS HELIOS 5 Diff.
23 02 03	Roche	COBAS HELIOS LMG
21 01 43	Roche	COBAS INTEGRA
23 02 02	Roche	COBAS MICROS
23 02 02	Roche	COBAS MINOS ST
23 02 02	Roche	COBAS MINOS STE
23 02 02	Roche	COBAS MINOS STEX
23 02 02	Roche	COBAS MINOS STX
23 02 02	Roche	COBAS MINOS VET
21 01 41	Roche	COBAS MIRA
21 01 41	Roche	COBAS MIRA PLUS
21 01 41	Roche	COBAS MIRA S
21 01 31	Roche	COBAS Ready (=SPOTCHEM)
23 02 04	Roche	COBAS VEGA
22 03 04	Abbott	Commander / FPC
21 02 01	AVL	Compact 1
21 02 01	AVL	Compact 2
21 03 03	Helena	Concept 2000
24 01 04	Coral Biotechnology	Coral Wave 120
21 01 19	Beckman	Creatinine-Analyzer 2
21 03 03	Eppendorf (JOKOO)	CTE 1000
21 03 03	Eppendorf (JOKOO)	CTE 5000
23 03 01	Ortho Diag.Syst. (J&J)	Cytoron Absolute
21 01 43	Bayer	DAX 24
21 01 43	Bayer	DAX 48
21 01 44	Bayer	DAX 72
21 01 44	Bayer	DAX 96
21 01 12	Bayer	DCA 2000
22 03 02	Wallac / Pharmacia	Delfia

EDMA CODE	Company	Instrument-Name
21 01 42	Kone	Delta
21 01 12	Bio-Rad	DIAMAT
22 03 04	Murex	DIAS System
21 01 41	Dade (Dupont)	Dimension 380
21 01 43	Dade (Dupont)	Dimension 760
21 01 43	Dade (Dupont)	Dimension AR
21 01 41	Dade (Dupont)	Dimension ES
21 01 43	Dade (Dupont)	Dimension XL
22 02 02	Bayer	DPA-1
21 01 01	Beckman	DU 5
21 01 42	Kone	Dynamic
21 01 52	Beckman	E2A
21 01 41	Eppendorf	EASY
21 01 31	Merck	EASY ST
21 01 11	Eppendorf	EBIO Compact
21 01 29	Eppendorf	EBIO Plus
22 03 01	Merck	Eclair
21 01 01	Merck	Eclipse
21 01 01	Merck	Eclipse Plus
21 01 01	Eppendorf	ECOM 6122
22 03 01	Behring	ECOM 6123 D
22 03 01	Eppendorf	ECOM 6123 D
21 01 01	Eppendorf	ECOM 6124 F
21 01 41	Eppendorf	ELAN
22 03 02	Boehringer Mannheim	ElecSys 1010
22 03 03	Boehringer Mannheim	ElecSys 2010
23 04 03	Ortho Diagnostic Systems (MLA)	Electra 1000 C
23 04 03	Ortho Diagnostic Systems (MLA)	Electra 1400 C
23 04 03	Ortho Diagnostic Systems (MLA)	Electra 1600 C
23 04 01	Ortho Diagnostic Systems (MLA)	Electra 700
23 04 02	Ortho Diagnostic Systems (MLA)	Electra 800
23 04 02	Ortho Diagnostic Systems (MLA)	Electra 900
23 04 02	Ortho Diagnostic Systems (MLA)	Electra 900 C
22 03 01	Behring	Elisa-Photometer
22 03 01	Behring	Elisa-Reader 301
22 03 01	Behring	Elisa-Reader 311
21 01 53	Beckman	ELISE
22 03 04	BIO-TEK	ELx808
21 01 32	Eppendorf	EPAC 6140
23 03 02	Coulter	Epics Elite
23 03 01	Coulter	Epics Profile II
21 01 41	Eppendorf	EPOS
21 01 42	Abbott	EPx
21 01 42	Eppendorf	ERIS
22 03 01	Boehringer Mannheim	ES 11
22 03 01	Boehringer Mannheim	ES 22
22 03 02	Boehringer Mannheim	ES 300
22 03 01	Boehringer Mannheim	ES 33
22 03 03	Boehringer Mannheim	ES 600
22 03 03	Boehringer Mannheim	ES 700
24 01 03	Difco	ESP
21 01 42	Eppendorf	ETEC
22 03 03	Sorin	ETI System
22 03 02	Behring	ETS Plus

EDMA CODE	Company	Instrument-Name
21 01 41	Chiron (Ciba Corning)	EXPRESS
23 02 01	Sysmex	F-520
23 02 01	Sysmex	F-820
23 03 01	Becton-Dickinson	FACSCalibur
23 03 01	Becton-Dickinson	FACScan
23 03 01	Becton-Dickinson	FACSCount
23 03 02	Becton-Dickinson	FACSort
23 03 02	Becton-Dickinson	FACStar ^{Plus}
23 03 01	Becton-Dickinson	FACStrak
23 03 02	Becton-Dickinson	FACSVantage
21 01 42	Menarini	Falcor 300
21 01 42	Menarini	Falcor 600
21 03 01	COSMO	FED 3200
22 02 01	Bio-Whittaker	Fiax
23 04 01	Behring	Fibrintimer 10
23 04 01	Behring	Fibrintimer II
23 04 01	Roche	Fibrometer
23 04 01	Labsystems	FP 900
23 04 01	Labsystems	FP 910
21 03 02	Olympus	Fractoscan
21 03 02	Olympus	Fractoscan Junior
21 02 02	Instrument. Lab. (IL)	GEM 6+
21 02 02	Instrument. Lab. (IL)	GEM PREMIER
21 02 02	Instrument. Lab. (IL)	GEM PREMIER PLUS
21 02 02	Instrument. Lab. (IL)	GEM STAT
26 02 02	Tecan	GENESIS RSP 200
23 02 02	SEAC	GENIUS
21 01 11	Beckman	Glucose-Analyzer 2
21 01 12	Chiron (Ciba Corning)	Glycomat
23 02 04	Bayer	H 1 E
23 02 04	Bayer	H 2
23 02 04	Bayer	H 3 RTC
23 02 04	Bayer	H 3 RTX
23 02 03	Bayer	H1
23 02 02	Bayer	H1JR
21 01 12	Menarini	HA 8110
21 01 12	Menarini	HA 8121
21 01 12	Menarini	HA 8140
23 04 03	BioMerieux	Hemolab
21 01 44	Boehringer Mannheim	HIT 737
21 01 44	Boehringer Mannheim	HIT 747 Series
21 01 12	Merck	Hitachi L 9000
21 01 41	Boehringer Mannheim	Hitachi 704
21 01 41	Boehringer Mannheim	Hitachi 705
21 01 43	Boehringer Mannheim	Hitachi 717
21 01 41	Boehringer Mannheim	Hitachi 904
21 01 42	Boehringer Mannheim	Hitachi 911
21 01 43	Boehringer Mannheim	Hitachi 917
21 03 03	Olympus	HITE SYSTEM 310
21 03 03	Olympus	HITE SYSTEM 320
21 03 03	Olympus	HITE SYSTEM 620
21 03 03	Menarini	HS 7
21 03 03	Sebia	HYDRASYS
22 02 01	Beckman	ICS

EDMA CODE	Company	Instrument-Name
21 01 01	Vitatron	IFP
21 02 01	Instrument. Lab. (IL)	IL 1304
21 02 01	Instrument. Lab. (IL)	IL 1306
21 02 01	Instrument. Lab. (IL)	IL 1312
21 02 02	Instrument. Lab. (IL)	IL 1400
21 02 02	Instrument. Lab. (IL)	IL 1405
21 02 01	Instrument. Lab. (IL)	IL 1420
21 02 02	Instrument. Lab. (IL)	IL 1430
21 02 01	Instrument. Lab. (IL)	IL 1610
21 02 01	Instrument. Lab. (IL)	IL 1620
21 02 02	Instrument. Lab. (IL)	IL 1630
21 02 02	Instrument. Lab. (IL)	IL 1640
21 02 02	Instrument. Lab. (IL)	IL 1650
21 02 07	Instrument. Lab. (IL)	IL 482
21 01 52	Instrument. Lab. (IL)	IL 614, 644, 654
21 02 07	Instrument. Lab. (IL)	IL 682
21 01 51	Instrument. Lab. (IL)	IL 943
21 02 01	Instrument. Lab. (IL)	IL BG3
21 01 44	Instrument. Lab. (IL)	ILab 1800
21 01 43	Instrument. Lab. (IL)	ILab 900
22 02 02	Beckman	IMMAGE
22 03 02	DPC	Immulite
22 03 03	Bayer	Immuno-I
22 03 02	Abbott	IMx / IMx Select
21 02 02	Diametrics	IRMA
21 01 01	Vitatron	ISP
23 02 03	Coulter	JR
23 02 03	Coulter	JS
23 02 02	Sysmex	K-1000
23 02 03	Sysmex	K-4500
23 02 02	Sysmex	K-800
23 04 01	Amelung	KC 10
23 04 01	Amelung	KC 10A
23 04 01	Amelung	KC 4
23 04 03	Amelung	KC 40
22 02 02	Boehringer Mannheim	KeySys
23 04 02	Ortho Diagnostic Systems (MLA)	Koagulab 16-605
23 04 01	Ortho Diagnostic Systems (MLA)	Koagulab CTS
21 01 31	Cholestech	L.D.X
22 03 03	Biochem Immunosystems	Labotec
21 01 01	Perkin-Elmer	Lambda 11
21 01 01	Perkin-Elmer	Lambda 12
21 01 01	Perkin-Elmer	Lambda 14
22 03 05	Abbott	LCx
22 03 02	Byk-Sangtek	LIA-mat S 300
22 03 01	Merck	LM-8
21 01 01	Lange (Dr. Lange)	LP 1
21 01 01	Lange (Dr. Lange)	LP 3
21 01 01	Lange (Dr. Lange)	LP 400
22 03 01	Sanofi Diagnostics Pasteur	LP 400
21 01 01	Lange (Dr. Lange)	LP 500
21 01 01	Lange (Dr. Lange)	LP 6A
23 02 01	Coulter	M430
23 02 01	Coulter	M530

EDMA CODE	Company	Instrument-Name
22 03 03	Merck (Biotrol)	MAGIA 7000
22 03 03	Merck (Biotrol)	MAGIA 8000
22 03 01	Chiron (Ciba Corning)	Magic Lite I
22 03 02	Chiron (Ciba Corning)	Magic Lite II
22 03 02	Chiron (Ciba Corning)	Magic Lite II+
22 03 02	Chiron (Ciba Corning)	Magic Lite SQ Allergy System
24 01 02	Mast	Mastascan
22 03 01	Abbott	Matrix
22 03 04	DPC / Biermann	MAX
23 02 04	Coulter	MAXM
23 04 02	Biodata	MCA 110
23 02 02	Coulter	MD 16
23 02 01	Coulter	MD 8
23 02 02	Coulter	MD II
21 01 43	Merck	MEGA
22 03 04	Murex	MIAS System
22 03 01	Murex	Micro ELISA MRX Reader
22 03 01	Murex	Micro ELISA Reader 7000
22 03 01	Organon Teknika	Microelisa Reader 100
22 03 01	Organon Teknika	Microelisa System
22 03 04	Hamilton	Microlab F.A.M.E.
24 01 02	Dade	Microscan Walkaway 40
24 01 02	Dade	Microscan Walkaway 96
22 03 01	Behring	Microtrack Reader SE System
21 05 03	Boehringer Mannheim	Miditron
21 05 03	Boehringer Mannheim	Miditron Junior
22 03 02	DPC	Milenia Max
22 03 02	BioMerieux	Mini-VIDAS
21 05 01	Menarini	MiniAution
21 01 01	Lange (Dr. Lange)	Miniphotometer 8
21 05 02	Boehringer Mannheim	Minitron
22 03 04	Merck	MIOS
21 01 01	Bayer	ML 1
21 01 01	Bayer	ML2
21 01 01	Bayer	ML3
21 01 01	Beckman	Model 42
21 01 41	Instrument. Lab. (IL)	Monarch 1000
21 01 42	Instrument. Lab. (IL)	Monarch 2000
21 01 42	Instrument. Lab. (IL)	Monarch PLUS
22 03 04	Organon Teknika	MPP 3000
22 03 05	Organon Teknika	NASBA QR System
23 02 03	Sysmex	NE-1500
23 02 04	Sysmex	NE-8000
21 03 01	Helena	NN
21 03 02	Helena	NN
21 03 01	Sebia	NN
21 03 02	Sebia	NN
21 05 04	Wako	NN
21 05 04	Yeong Deong	NN
21 01 42	Nova	Nucleus
21 02 04	AVL	OMNI
22 03 04	Sorin	Omni System
23 02 03	Coulter	Onyx
21 01 42	Bayer	OPERA

EDMA CODE	Company	Instrument-Name
21 02 01	AVL	OPTI 1
21 01 42	Coulter	Optichem 500
21 01 42	Kone	Optima
23 04 01	BioMerieux	Option 2-8
22 03 02	Behring	OPUS
22 03 03	Behring	OPUS Magnum
22 03 02	Behring	OPUS PLUS
21 02 07	Radiometer	OSM 3
23 04 04	Biodata	PAP 4
21 03 01	Beckman	PARAGON
21 03 03	Beckman	PARAGON CZE 2000
21 01 43	Dade (Dupont)	Paramax-RX
23 03 02	Partec	PAS-II
24 01 02	Difco	Pasco
21 01 01	Eppendorf	PCP 6121
21 01 43	American Monitor	Perspective/Excel
23 04 04	Ortho Diagnostic Systems / Dade	PFA-100
21 01 01	Boehringer Mannheim	Photometer 4010
21 01 01	Boehringer Mannheim	Photometer 4020
22 03 02	Hybritech	Photon Era
22 03 01	Hybritech	Photon II
23 05 01	Olympus	PK 7100
23 05 01	Olympus	PK 7200
23 05 01	Olympus	PK 80
22 03 03	Rosys	Plato 1300
22 03 04	Rosys	Plato 3300
21 01 01	Riele	PM 210
21 01 01	Riele	PM 310
21 01 01	Riele	PMC 2000
26 02 02	Abbott	PPC
26 02 09	Boehringer Mannheim	PPV
21 03 03	Helena	PRINCE-CE
22 03 04	Abbott	PRISM
22 03 05	Becton Dickinson	ProbeTec
21 01 41	Kone	Progress
21 01 41	Kone	Progress Plus
21 01 01	Unicam	PU 8610
21 01 01	Unicam	PU 8620
23 02 01	Becton-Dickinson	QBC Autoread
23 02 01	Becton-Dickinson	QBC I
23 02 01	Becton-Dickinson	QBC II
22 02 02	Sanofi Diagnostics Pasteur	QM 300
22 03 01	Abbott	Quantum
22 03 01	Abbott	Quantum II
22 03 03	Kone	Quartus
23 05 01	Lifescreeen	QUASAR 240
23 05 01	Lifescreeen	QUASAR 400
24 01 04	Difco	Questor
22 02 01	Orion Diagnostica	QuikRead 3
23 03 01	Sysmex	R-2000/3000
23 02 11	Sysmex	R-3000
21 01 41	Bayer	RA 1000
21 01 41	Bayer	RA 50
21 01 41	Bayer	RA 500

EDMA CODE	Company	Instrument-Name
21 01 41	Bayer	RA XT
21 01 42	Bayer	RA2000
22 03 04	Bio-Rad	Radius
24 01 02	Mast	Radius
24 01 04	Orbec	Ramus
21 05 03	Behring	RAPIMAT
21 01 31	Boehringer Mannheim	Reflotron
21 01 31	Boehringer Mannheim	Reflotron Sprint
21 03 03	Helena	REP
23 02 02	Coulter	S770
23 02 03	Coulter	S880
24 01 02	Becton Dickinson	Sceptor
24 01 02	Becton Dickinson	Sceptorpette
23 02 04	Sysmex	SE-9000
21 01 41	Merck	Selectra
24 01 02	Radiometer	Sensititre
24 01 03	Difco	Sentinel
21 01 31	Bayer	Seralyzer
22 03 01	Biochem Immunosystems	Serazyme I
22 03 01	Biochem Immunosystems	Serazyme II
23 02 04	Sysmex	SF-3000
22 03 02	Behring	Solaris
22 03 04	Sanofi Diagnostics Pasteur	Sonata
21 01 41	Kone	Specific
21 01 41	Kone	Specific Supra
21 01 01	Milton-Roy	Spectronic 1001+
21 01 01	Milton-Roy	Spectronic 1101
21 01 01	Milton-Roy	Spectronic 1201
21 01 01	Milton-Roy	Spectronic 20 / 20 D
21 01 01	Milton-Roy	Spectronic 21 DV
21 01 01	Milton-Roy	Spectronic 21 MV
21 01 01	Milton-Roy	Spectronic 21 UVD
21 01 01	Milton-Roy	Spectronic 3000
21 01 01	Milton-Roy	Spectronic 301
21 01 01	Milton-Roy	Spectronic 401
21 01 01	Milton-Roy	Spectronic 501
21 01 01	Milton-Roy	Spectronic 601
21 01 41	Abbott	Spectrum
21 01 31	Eppendorf	SPOTCHEM
21 01 31	Menarini	SPOTCHEM
22 03 02	Biochem Immunosystems	SR 1
23 02 04	Coulter	ST
23 04 03	Boehringer Mannheim / Stago	STA
23 04 01	Stago	STA 2
23 04 01	Stago	STA 4
23 04 01	Stago	STA 8
23 04 01	Stago	STA 888
23 04 03	Boehringer Mannheim / Stago	STA Compact
21 02 01	Nova	Stat Profile 3
21 02 02	Nova Biomedical	Stat Profile Plus 9,10
21 02 02	Nova Biomedical	Stat Profile1,2,4,5,7
23 02 04	Coulter	STKR
23 02 04	Coulter	STKS
23 04 03	Stratec	Stracomat

EDMA CODE	Company	Instrument-Name
22 03 02	Dade	Stratus I
22 03 02	Dade	Stratus II
24 01 02	I2R	STRSCAN
22 03 04	Ortho Diagnostic Systems	Summit
26 02 02	Ortho Diag.Syst. (J&J)	Summit Plus 2
26 02 02	Ortho Diag.Syst. (J&J)	Summit Processor
21 05 04	Menarini	SuperAution
21 05 04	Boehringer Mannheim	Supertron
21 01 42	Beckman	Synchron CX 3 DELTA
21 01 42	Beckman	Synchron CX 4
21 01 42	Beckman	Synchron CX 4 DELTA
21 01 42	Beckman	Synchron CX 5
21 01 42	Beckman	Synchron CX 5 DELTA
21 01 43	Beckman	Synchron CX 7
21 01 43	Beckman	Synchron CX 7 DELTA
21 01 41	Beckman	Synchron CX3
21 02 04	Instrument. Lab. (IL)	Synthesis 25 /35
22 03 01	Dynatech	System 5000
22 03 04	Dynatech	System 7000
23 02 02	Coulter	T540
23 02 02	Coulter	T660
23 02 03	Coulter	T890
22 03 02	Abbott	TDx / TDxFLx
22 03 03	Organon Teknika	Tektime
23 04 03	Behnk	Thrombolysar
23 04 01	Nycomed	Thrombotrack
24 01 02	Dade	Touchscan
22 02 01	Behring	TurbiTimer
22 02 01	Orion Diagnostica	Turbox
24 01 04	Sysmex	UF-100
21 01 42	Kone	Ultra
22 03 02	Pharmacia	UniCAP
24 01 02	BioMerieux	Uniscept
21 01 19	Beckman	Urea-Analyzer 2
21 05 02	Boehringer Mannheim	Urilux
21 05 02	IMACO	URIMAC
24 01 04	BioMerieux	URIMAT
21 05 03	Boehringer Mannheim	Urotron RL 9
21 01 12	Bio-Rad	VARIANT
22 03 02	BioMerieux	VIDAS
24 01 02	Biokit	Videobact
21 01 31	Abbott	Vision
24 01 03	BioMerieux	Vital 1200
24 01 03	BioMerieux	Vital 200
24 01 02	BioMerieux	Vitek 120
24 01 02	BioMerieux	Vitek 240
24 01 02	BioMerieux	Vitek 32
24 01 02	BioMerieux	Vitek 60
24 01 02	BioMerieux	Vitek Junior
22 03 03	Ortho Diag.Syst. (J&J)	VITROS ECI
21 01 42	Ortho Diag.Syst. (J&J)	VITROS(Ektachem) 250
21 01 42	Ortho Diag.Syst. (J&J)	VITROS(Ektachem) 500
21 01 43	Ortho Diag.Syst. (J&J)	VITROS(Ektachem) 750 XRC
21 01 43	Ortho Diag.Syst. (J&J)	VITROS(Ektachem) 950 IRC

EDMA CODE	Company	Instrument-Name
21 01 31	Ortho Diag.Syst. (J&J)	VITROS(Ektachem) DT 60
21 01 41	Abbott	VP
21 01 11	Yellow Springs	YSI
23 02 01	Coulter	ZM